

---

## APC Anti-mouse CD150 (TC15-12F12.2)

### Product Information

Catalog No.:	APC-09790
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	TC15-12F12.2
IsoType:	IgG2a, lambda
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	CDw150, Slamf1, Slam, ESTM51, IPO-3, 4933415F16, AA177906
Uniprot ID :	Q9QUM4

### Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution: Volume per test:	5 µ L