

## FITC Anti-Human CD81 (1.3.3.22)

### Product Information

Catalog No.:	FITC-30141
Conjugation:	FITC
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1.3.3.22
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	TSPAN28 26 kDa cell surface protein TAPA-1 APA1 CD81 CD81 antigen Target of the antiproliferative antibody 1 Tetraspanin-28 Tspan-28
Uniprot ID :	P60033

### Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: