

---

## AP-conjugated Affinipure Goat Anti- Rat IgG (H+L)

### Product Information

Catalog No.:	FNSA-0028
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	Alkaline Phosphatase purified from calf intestine
Stabilizer:	5mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application: ELISA, WB

Suggested dilution: 1:1000-1:3000 for ELISA and Western blotting

**Note: This product is for research use only, not for diagnostic or therapeutic use.**