

---

## FITC Conjugated Donkey Anti Rabbit IgG (H+L)

### Product Information

Catalog No.:	FNSA-0034
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was affinity isolated.
Fluorophore:	Fluorescein-5-isothiocyanate(FITC) Amax=492nm;Emax=520nm
Fluorophore/Protein:	10.4ug/mg; 3.9 moles FITC per mole IgG
Stabilizer:	5mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application: IF, Flow cyt  
Suggested dilution: 1:20-1:100 for most applications

**Note: This product is for research use only, not for diagnostic or therapeutic use.**