

---

## Cy3–Rabbit Anti-Goat IgG (H+L)

### Product Information

Catalog No.:	FNSA-0080
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	Purified from seaweed
Fluorophore:	Cy3 bis-NHS ester $A_{max}=550nm, E_{max}=570nm$
Fluorophore/Protein Absorbance Ratio:	$A_{552}/A_{280}=2.2$
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application:	IF, FC
Recommended dilution:	1:20-1:100

**Note: This product is for research use only, not for diagnostic or therapeutic use.**