
HRP conjugated Goat anti Monkey IgM

Product Information

| | |
|--------------|--|
| Catalog No.: | FNSA-0100 |
| Size: | 100µl/500µl |
| Form: | liquid |
| Conjugation: | HRP purified from horseradish root |
| Purity: | The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Stabilizer: | 3mg/ml BSA (IgG-Free, Protease-Free) |
| Storage: | 0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20°C for 24 months , avoid repeated freezing and thawing. |

Application

Tested Application: ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.