
Poly-HRP conjugated Mouse Anti-human IgG

Product Information

Catalog No.:	FNSA-0114
Size:	100µL/500µL
Form:	liquid
Conjugation:	HRP purified from horseradish root
Stabilizer:	3mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, and 50% glycerol, pH7.6. -20°C for 24 months.

Application

Tested Application:ELISA

Suggested dilution: 1:5000-1:20000

Note: This product is for research use only, not for diagnostic or therapeutic use.