

---

## HRP-Goat Anti-Mouse IgG ( $\gamma$ -chain specific)

### Product Information

Catalog No.:	FNSA-0161
Size:	100 $\mu$ l/500 $\mu$ l
Form:	liquid
Conjugation:	HRP purified from horseradish root
Immunogen:	Prepared from purified mouse IgG $\gamma$ chain
Stabilizer:	3mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application:ELISA

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

**Note: This product is for research use only, not for diagnostic or therapeutic use.**