

AOX1 antibody

Product Information

Catalog No.:	FNab00281
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	aldehyde oxidase 1
Synonyms:	Aldehyde oxidase Aldehyde oxidase 1 Azaheterocycle hydroxylase AOX1 AO
Observed MW:	148 kDa
Uniprot ID :	Q06278

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IP: 1:200-1:1000; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab00281(AOX1 antibody) at dilution of 1:100



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab00281(AOX1 antibody) at dilution of 1:300