

## ANKRD33 antibody

### Product Information

Catalog No.:	FNab00412
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	ankyrin repeat domain 33
Synonyms:	Photoreceptor ankyrin repeat protein Ankyrin repeat domain-containing protein 33 ANKRD33 C12orf7 PANKY
Observed MW:	49 kDa
Uniprot ID :	Q7Z3H0

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00412(ANKRD33 Antibody) at dilution of 1:50



Y79 cells were subjected to SDS PAGE followed by western blot with FNab00412(ANKRD33 Antibody) at dilution of 1:600