

BCAT1 antibody

Product Information

Catalog No.:	FNab00830
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	branched chain aminotransferase 1, cytosolic
Synonyms:	Branched-chain-amino-acid aminotransferase, cytosolic (BCAT(c)) Protein ECA39 BCAT1 BCT1 ECA39
Observed MW:	43-45 kDa
Uniprot ID :	P54687

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using FNab00830(BCAT1 Antibody) at dilution of 1:50



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab00830(BCAT1 antibody) at dilution of 1:300