
BIN3 antibody

Product Information

Catalog No.:	FNab00900
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in cytokinesis and septation where it has a role in the localization of F-actin.

Immunogen information

Immunogen:	bridging integrator 3
Synonyms:	Bridging integrator 3 BIN3
Observed MW:	30 kDa
Uniprot ID :	Q9NQY0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab00900(BIN3 antibody) at dilution of 1:300