

BCAN antibody

Product Information

Catalog No.:	FNab00954
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	brevican
Synonyms:	Brevican core protein Brain-enriched hyaluronan-binding protein (BEHAB) Chondroitin sulfate proteoglycan 7 BCAN BEHAB CSPG7
Observed MW:	145 kDa, 80 kDa
Uniprot ID :	Q96GW7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human cerebellum tissue slide using FNab00954(Brevican antibody at dilution of 1:200



human brain tissue were subjected to SDS PAGE followed by western blot with FNab00954(Brevican antibody) at dilution of 1:500