

NRDE2 antibody

Product Information

Catalog No.:	FNab01018
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	chromosome 14 open reading frame 102
Synonyms:	Nuclear exosome regulator NRDE2 Protein NRDE2 homolog NRDE2 C14orf102
Observed MW:	133 kDa
Uniprot ID :	Q9H7Z3

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	



IP Result of anti-C14orf102 (IP:FNab01018, 5ug;
Detection:FNab01018 1:500) with MCF-7 cells
lysate 1200ug.



MCF-7 cells were subjected to SDS PAGE
followed by western blot with
FNab01018(C14orf102 Antibody) at dilution of
1:600