

CEP57L1 antibody

Product Information

Catalog No.:	FNab01128
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	chromosome 6 open reading frame 182
Synonyms:	Centrosomal protein CEP57L1 Centrosomal protein 57kDa-like protein 1 Centrosomal protein of 57 kDa-related protein (Cep57R) Cep57-related protein CEP57L1 C6orf182 CEP57R
Observed MW:	55-65 kDa
Uniprot ID :	Q8IYX8

Application

Reactivity:	Human, Mouse, Monkey
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01128(C6orf182 Antibody) at dilution of 1:600