

## PLGRKT antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab01145  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information

|              |  |
|--------------|--|
| Immunogen:   | chromosome 9 open reading frame 46                   |
| Synonyms:    | Plasminogen receptor (KT) (Plg-R(KT)) PLGRKT C9orf46 |
| Observed MW: | 17-20 kDa  |
| Uniprot ID : | Q9HBL7   |

### Application

|                           |              |
|---------------------------|--------------|
| Reactivity:               | Human        |
| Tested Application:       | ELISA, WB    |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image:                    |              |



A431 cells were subjected to SDS PAGE followed by western blot with FNab01145( C9orf46 Antibody) at dilution of 1:1000