

## CADM4 antibody

### Product Information

Catalog No.:	FNab01192
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the cell-cell adhesion. Has calcium-and magnesium-independent cell-cell adhesion activity. May have tumor-suppressor activity.

### Immunogen information

Immunogen:	cell adhesion molecule 4
Synonyms:	Cell adhesion molecule 4 Immunoglobulin superfamily member 4C (IgSF4C) Nectin-like protein 4 (NECL-4) TSLC1-like protein 2 CADM4 IGSF4C NECL4 TSL2
Observed MW:	65-70 kDa
Uniprot ID :	Q8NFZ8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab01192(CADM4 antibody) at dilution of 1:1000