

CAMLG antibody

Product Information

Catalog No.:	FNab01240
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Likely involved in the mobilization of calcium as a result of the TCR/CD3 complex interaction.
Binds to cyclophilin B.

Immunogen information

Immunogen:	calcium modulating ligand
Synonyms:	Guided entry of tail-anchored proteins factor CAMLG Calcium signal-modulating cyclophilin ligand CAMLG CAML GET2
Observed MW:	36 kDa
Uniprot ID :	P49069

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01240(CAMLG antibody) at dilution of 1:1000