

CTSD antibody

Product Information

Catalog No.:	FNab01306
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	cathepsin D
Synonyms:	Cathepsin D Cathepsin D light chain Cathepsin D heavy chain CTSD CPSD
Observed MW:	48-52 kDa
Uniprot ID :	P07339

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab01306 (Cathepsin D Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab01306 (CTSD antibody) at dilution of 1:1000