

## CCDC42 antibody

### Product Information

Catalog No.:	FNab01357
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	coiled-coil domain containing 42
Synonyms:	Coiled-coil domain-containing protein 42 CCDC42 CCDC42A
Observed MW:	38 kDa
Uniprot ID :	Q96M95

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab01357(CCDC42 antibody) at dilution of 1:100