

## KLRD1 antibody

### Product Information

Catalog No.:	FNab01509
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	killer cell lectin-like receptor subfamily D, member 1
Synonyms:	Natural killer cells antigen CD94 KP43 Killer cell lectin-like receptor subfamily D member 1 NK cell receptor KLRD1 CD94
Observed MW:	25-30 kDa
Uniprot ID :	Q13241

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01509(KLRD1 Antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01509(KLRD1 antibody) at dilution of 1:500