

## CDK9 antibody

### Product Information


Catalog No.:	FNab01564
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	cyclin-dependent kinase 9
Synonyms:	Cyclin-dependent kinase 9 C-2K Cell division cycle 2-like protein kinase 4 Cell division protein kinase 9 Serine/threonine-protein kinase PITALRE Tat-associated kinase complex catalytic subunit CDK9 CDC2L4 TAK
Observed MW:	43 kDa
Uniprot ID :	P50750

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab01564(CDK9 antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab01564(CDK9 antibody) at dilution of 1:1000