

SH3KBP1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab01712 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|---|
| Immunogen: | SH3-domain kinase binding protein 1 |
| Synonyms: | SH3 domain-containing kinase-binding protein 1 CD2-binding protein 3 (CD2BP3) Cbl-interacting protein of 85 kDa Human Src family kinase-binding protein 1 (HSB-1) SH3KBP1 CIN85 |
| Observed MW: | 85 kDa |
| Uniprot ID : | Q96B97 |

Application

| | |
|---------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution: WB: | 1:500-1:2000 |
| Image: | |



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab01712(CIN85 Antibody) at dilution of 1:1000