

CLEC10A antibody

Product Information

Catalog No.:	FNab01745
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	C-type lectin domain family 10, member A
Synonyms:	C-type lectin domain family 10 member A C-type lectin superfamily member 14 Macrophage lectin 2 CLEC10A CLECSF13 CLECSF14 HML
Observed MW:	32 kDa
Uniprot ID :	Q8IUN9

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab01745(CLEC10A/CD301 antibody) at dilution of 1:500