

CLIP4 antibody

Product Information

Catalog No.:	FNab01765
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	CAP-GLY domain containing linker protein family, member 4
Synonyms:	CAP-Gly domain-containing linker protein 4 Restin-like protein 2 CLIP4 RSNL2
Observed MW:	65-75 kDa
Uniprot ID :	Q8N3C7

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Y79 cells were subjected to SDS PAGE followed by western blot with FNab01765 (CLIP4 antibody) at dilution of 1:300