

MYB antibody

Product Information


Catalog No.:	FNab01789
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Transcriptional activator Myb
Synonyms:	MYB Proto-oncogene c-Myb c-Myb cMyb Transcriptional activator Myb
Observed MW:	72-80 kDa
Uniprot ID :	P10242

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab01789(c-Myb antibody) at dilution of 1:500