

## COG1 antibody

### Product Information

Catalog No.:	FNab01828
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Required for normal Golgi function.

### Immunogen information


Immunogen:	component of oligomeric golgi complex 1
Synonyms:	Conserved oligomeric Golgi complex subunit 1 (COG complex subunit 1) Component of oligomeric Golgi complex 1 COG1 KIAA1381 LDLB
Observed MW:	109 kDa
Uniprot ID :	Q8WTW3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab01828(COG1 antibody) at dilution of 1:100



U251 cell lysates were subjected to SDS PAGE followed by western blot with FNab01828(COG1 antibody) at dilution of 1:1000