

## CFH antibody

### Product Information

Catalog No.:	FNab01855
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Factor H functions as a cofactor in the inactivation of C3b by factor I and also increases the rate of dissociation of the C3bBb complex(C3 convertase) and the(C3b)NBB complex(C5 convertase) in the alternative complement pathway.

### Immunogen information


Immunogen:	complement factor H
Synonyms:	Complement factor H H factor 1 CFH HF HF1 HF2
Observed MW:	150 kDa
Uniprot ID :	P08603

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab01855(CFH Antibody) at dilution of 1:50



human blood tissue were subjected to SDS PAGE followed by western blot with FNab01855(CFH Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours