

CPEB4 antibody

Product Information

Catalog No.:	FNab01918
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cytoplasmic polyadenylation element binding protein 4
Synonyms:	Cytoplasmic polyadenylation element-binding protein 4 (CPE-BP4, CPE-binding protein 4, hCPEB-4) CPEB4 KIAA1673
Observed MW:	80 kDa
Uniprot ID :	Q17RY0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab01918(CPEB4 Antibody) at dilution of 1:600