

CWC25 antibody

Product Information

Catalog No.:	FNab02086
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	coiled-coil domain containing 49
Synonyms:	Pre-mRNA-splicing factor CWC25 homolog Coiled-coil domain-containing protein 49 Spliceosome-associated protein homolog CWC25 CWC25 CCDC49
Observed MW:	50 kDa
Uniprot ID :	Q9NXE8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab02086(CCDC49 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02086(CCDC49 Antibody) at dilution of 1:1000