

CYP11A1 antibody

Product Information

Catalog No.:	FNab02143
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the side-chain cleavage reaction of cholesterol to pregnenolone.

Immunogen information


Immunogen:	cytochrome P450, family 11, subfamily A, polypeptide 1
Synonyms:	Cholesterol side-chain cleavage enzyme, mitochondrial CYPXIA1 Cholesterol desmolase Cytochrome P450 11A1 Cytochrome P450(scc) CYP11A1 CYP11A
Observed MW:	49-52 kDa
Uniprot ID :	P05108

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab02143(CYP11A1 Antibody) at dilution of 1:200



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab02143(CYP11A1 antibody) at dilution of 1:800