

NR0B1 antibody

Product Information

Catalog No.:	FNab02248
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nuclear receptor subfamily 0, group B, member 1
Synonyms:	Nuclear receptor subfamily 0 group B member 1 DSS-AHC critical region on the X chromosome protein 1 Nuclear receptor DAX-1 NR0B1 AHC DAX1
Observed MW:	54 kDa
Uniprot ID :	P51843

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab02248(DAX-1 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab02248(DAX-1 Antibody) at dilution of 1:300