

DDA1 antibody

Product Information

Catalog No.:	FNab02282
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins.

Immunogen information

Immunogen:	DET1 and DDB1 associated 1
Synonyms:	DET1- and DDB1-associated protein 1 Placenta cross-immune reaction antigen 1 (PCIA-1) DDA1 C19orf58 PCIA1
Observed MW:	12 kDa
Uniprot ID :	Q9BW61

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab02282(DDA1 antibody) at dilution of 1:100