

DHX16 antibody

Product Information

Catalog No.:	FNab02376
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable ATP-binding RNA helicase involved in pre-mRNA splicing.

Immunogen information


Immunogen:	DEAH(Asp-Glu-Ala-His) box polypeptide 16
Synonyms:	Pre-mRNA-splicing factor ATP-dependent RNA helicase DHX16 ATP-dependent RNA helicase #3 DEAH-box protein 16 DHX16 DBP2 DDX16 KIAA0577 PRP2
Observed MW:	119 kDa
Uniprot ID :	O60231

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab02376(DHX16 antibody) at dilution of 1:10



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab02376(DHX16 antibody) at dilution of 1:1000