

DHX35 antibody

Product Information

Catalog No.:	FNab02379
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in pre-mRNA splicing.

Immunogen information


Immunogen:	DEAH(Asp-Glu-Ala-His) box polypeptide 35
Synonyms:	Probable ATP-dependent RNA helicase DHX35 DEAH box protein 35 DHX35 C20orf15 DDX35
Observed MW:	70 kDa
Uniprot ID :	Q9H5Z1

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab02379(DHX35 Antibody) at dilution of 1:50



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab02379(DHX35 Antibody) at dilution of 1:300