

DMWD antibody

Product Information

Catalog No.:	FNab02428
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Could have a regulatory function in meiosis.

Immunogen information

Immunogen:	dystrophia myotonica, WD repeat containing
Synonyms:	Dystrophia myotonica WD repeat-containing protein Dystrophia myotonica-containing WD repeat motif protein Protein 59 Protein DMR-N9 DMWD DM9
Observed MW:	70 kDa
Uniprot ID :	Q09019

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab02428(DMWD antibody) at dilution of 1:600