

DNAJC9 antibody

Product Information

Catalog No.:	FNab02470
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8.

Immunogen information

Immunogen:	DnaJ(Hsp40) homolog, subfamily C, member 9
Synonyms:	DnaJ homolog subfamily C member 9 (HDJC9) DnaJ protein SB73 DNAJC9
Observed MW:	30-35 kDa
Uniprot ID :	Q8WXX5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02470(DNAJC9 Antibody) at dilution of 1:1000