

## DPEP1 antibody

### Product Information

Catalog No.:	FNab02513
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Hydrolyzes a wide range of dipeptides. Implicated in the renal metabolism of glutathione and its conjugates. Converts leukotriene D4 to leukotriene E4; it may play an important role in the regulation of leukotriene activity.

### Immunogen information


Immunogen:	dipeptidase 1(renal)
Synonyms:	Dipeptidase 1 Beta-lactamase Dehydropeptidase-I Microsomal dipeptidase Renal dipeptidase (hRDP) DPEP1 MDP RDP
Observed MW:	46 kDa
Uniprot ID :	P16444

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab02513(DPEP1 Antibody) at dilution of 1:200



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab02513(DPEP1 Antibody) at dilution of 1:300