

DSN1 antibody

Product Information

Catalog No.:	FNab02546
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore formation during mitosis.

Immunogen information

Immunogen:	DSN1, MIND kinetochore complex component, homolog
Synonyms:	Kinetochore-associated protein DSN1 homolog DSN1 C20orf172 MIS13
Observed MW:	40 kDa
Uniprot ID :	Q9H410

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02546(DSN1 Antibody) at dilution of 1:1500