

DUS3L antibody

Product Information

Catalog No.:	FNab02560
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the synthesis of dihydrouridine, a modified base found in the D-loop of most tRNAs.

Immunogen information


Immunogen:	dihydrouridine synthase 3-like
Synonyms:	tRNA-dihydrouridine(47) synthase [NAD(P)(+)]-like mRNA-dihydrouridine synthase DUS3L tRNA-dihydrouridine synthase 3-like DUS3L
Observed MW:	72 kDa
Uniprot ID :	Q96G46

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab02560(DUS3L antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab02560(DUS3L antibody) at dilution of 1:500