

DUSP22 antibody

Product Information

Catalog No.:	FNab02568
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Activates the Jnk signaling pathway. Dephosphorylates and deactivates p38 and stress-activated protein kinase/c-Jun N-terminal kinase(SAPK/JNK)(By similarity).

Immunogen information


Immunogen:	dual specificity phosphatase 22
Synonyms:	Dual specificity protein phosphatase 22 JNK-stimulatory phosphatase-1 (JSP-1) Low molecular weight dual specificity phosphatase 2 (LMW-DSP2) Mitogen-activated protein kinase phosphatase x (MAP kinase phosphatase x, MKP-x) DUSP22 JSP1 LMWDSP2 MKPX
Observed MW:	21 kDa
Uniprot ID :	Q9NRW4

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab02568(DUSP22 antibody) at dilution of 1:50



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02568(DUSP22 Antibody) at dilution of 1:1000