

DYNC1LI2 antibody

Product Information

Catalog No.:	FNab02586
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in binding dynein to membranous organelles or chromosomes.


Immunogen information

Immunogen:	dynein, cytoplasmic 1, light intermediate chain 2
Synonyms:	Cytoplasmic dynein 1 light intermediate chain 2 Dynein light intermediate chain 2, cytosolic (LIC-2) LIC53/55 DYNC1LI2 DNCLI2 LIC2
Observed MW:	60-65 kDa
Uniprot ID :	O43237

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000

Image:



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab02586(DYNC1LI2 antibody) at dilution of 1:600