

## DYNLRB1 antibody

### Product Information

Catalog No.:	FNab02589
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules.


### Immunogen information

Immunogen:	dynein, light chain, roadblock-type 1
Synonyms:	Dynein light chain roadblock-type 1 Bithoraxoid-like protein (BLP) Dynein light chain 2A, cytoplasmic Dynein-associated protein Km23 Roadblock domain-containing protein 1 DYNLRB1 BITH DNCL2A DNLC2A ROBLD1
Observed MW:	20 kDa
Uniprot ID :	Q9NP97

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000

Image:



293T cells were subjected to SDS PAGE followed by western blot with FNab02589(DYNLRB1 antibody) at dilution of 1:1000