

## DZIP1 antibody

### Product Information

Catalog No.:	FNab02598
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


FNab02737

### Immunogen information


Immunogen:	DAZ interacting protein 1
Synonyms:	Cilium assembly protein DZIP1 DAZ-interacting protein 1/2 DAZ-interacting zinc finger protein 1 DZIP1 DZIP DZIP2 KIAA0996
Observed MW:	100-113 kDa
Uniprot ID :	Q86YF9

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-DZIP1 (IP:FNab02598, 4ug;  
Detection:FNab02598 1:300) with HEK-293 cells  
lysate 880ug.



HEK-293 cells were subjected to SDS PAGE  
followed by western blot with FNab02598(DZIP1  
antibody) at dilution of 1:300