

## EAPP antibody

### Product Information

Catalog No.:	FNab02610
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play an important role in the fine-tuning of both major E2F1 activities, the regulation of the cell-cycle and the induction of apoptosis. Promotes S-phase entry, and inhibits p14(ARP) expression.

### Immunogen information

Immunogen:	E2F-associated phosphoprotein
Synonyms:	E2F-associated phosphoprotein (EAPP) EAPP C14orf11
Observed MW:	40-45 kDa
Uniprot ID :	Q56P03

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab02610( EAPP Antibody) at dilution of 1:600