

MAPRE1 antibody

Product Information

Catalog No.:	FNab02611
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	microtubule-associated protein, RP/EB family, member 1
Synonyms:	Microtubule-associated protein RP/EB family member 1 APC-binding protein EB1 End-binding protein 1 (EB1) MAPRE1
Observed MW:	35 kDa
Uniprot ID :	Q15691

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab02611(MAPRE1 antibody) at dilution of 1:100



IP Result of anti-EB1 (IP:FNab02611, 4ug; Detection:FNab02611 1:1000) with HeLa cells lysate 1200ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02611(MAPRE1 antibody) at dilution of 1:500