

## EIF4B antibody

### Product Information

Catalog No.:	FNab02720
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


EIF4B (eukaryotic translation initiation factor 4B) is a protein-coding gene. GO annotations related to this gene include RNA binding and translation initiation factor activity. An important paralog of this gene is EIF4H.

### Immunogen information


Immunogen:	eukaryotic translation initiation factor 4B
Synonyms:	Eukaryotic translation initiation factor 4B (eIF-4B) EIF4B
Observed MW:	85 kDa
Uniprot ID :	P23588

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500 - 1:2000; IF: 1:50 - 1:200
Image:	



Immunofluorescence analysis of U2OS cells using FNab02720(EIF4B antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab02720(EIF4B antibody) at dilution of 1:1000