

ELP4 antibody

Product Information

Catalog No.:	FNab02749
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as subunit of the RNA polymerase II elongator complex, which is a histone acetyltransferase component of the RNA polymerase II(Pol II) holoenzyme and is involved in transcriptional elongation. Elongator may play a role in chromatin remodeling and is involved in acetylation of histones H3 and probably H4.

Immunogen information


Immunogen:	elongation protein 4 homolog
Synonyms:	Elongator complex protein 4 (hELP4) PAX6 neighbor gene protein ELP4 C11orf19 PAXNEB
Observed MW:	46 kDa
Uniprot ID :	Q96EB1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab02749(ELP4 antibody) at dilution of 1:50



human cerebellum tissue were subjected to SDS PAGE followed by western blot with FNab02749(ELP4 antibody) at dilution of 1:300