

EPB41L1 antibody

Product Information

Catalog No.:	FNab02796
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function to confer stability and plasticity to neuronal membrane via multiple interactions, including the spectrin-actin-based cytoskeleton, integral membrane channels and membrane-associated guanylate kinases.

Immunogen information

Immunogen:	erythrocyte membrane protein band 4.1-like 1
Synonyms:	Band 4.1-like protein 1 Erythrocyte membrane protein band 4.1-like 1 Neuronal protein 4.1 (4.1N) EPB41L1 KIAA0338
Observed MW:	100-130 kDa
Uniprot ID :	Q9H4G0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab02796(EPB41L1 Antibody) at dilution of 1:600