

LIMA1 antibody

Product Information

Catalog No.:	FNab02813
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1D0
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	LIM domain and actin binding 1
Synonyms:	LIM domain and actin-binding protein 1 Epithelial protein lost in neoplasm LIMA1 EPLIN SREBP3
Observed MW:	90 kDa, 110 kDa
Uniprot ID :	Q9UHB6

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab02813(LIMA1 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab02813(LIMA1 antibody) at dilution of 1:1000